

## Experiment 5: Microscopy: *Candida* morphological alterations

### Protocol 5.1: Set up cultures for *Candida* microscopy (Week 8)

1. You will be provided with a 10 mL *Candida albicans* culture in PIM medium. You will use this culture to set up 4 cultures with different concentrations of Antibiotic X (0, 10, 100, and 300  $\mu$ M X). Label 4 sterile universal tubes appropriately.
2. Aseptically transfer 1 mL of culture into the four labelled sterile universal tubes
3. Add 9 mL of fresh PIM medium to each tube.
4. Add the appropriate volumes of a 0.1 M Antibiotic X solution, so that you achieve final concentrations of 10, 100 and 300  $\mu$ M of X in three of the tubes. Leave the fourth tube without X as a negative control. Let the cultures grow overnight at 37°C.

**Table 5.1 Antibiotic X dilutions (prepared using 0.1 M antibiotic X)**

Tube	Final [Antibiotic X] ( $\mu$ M)	Amount of X to add
1	0	
2	10	
3	100	
4	300	

### Protocol 5.2: *Candida* microscopy (Week 10)

5. **In Week 10**, observe samples of each of your four cultures under the microscope. Were there any morphological alterations? Do you observe any other changes? Could you explain these changes?
6. Use a haemocytometer to count yeast and filamentous cells. Calculate the % of each morphological type. Record your results in Table 5.2.

**Table 5.2: Morphological alterations of *Candida albicans* after Antibiotic X treatment.**

Antibiotic concentration	X	Number of yeast cells	Number of hyphal cells